



# **MEDICAL CARE**

TEL: 0086-510-82123587

WEB: <https://www.aupcon.com>

ADD: No.5 Minshan Road Xinwu district, Wuxi China

**Aupcon Medical Technology (wuxi) Co., Ltd**

# COMPANY PROFILE

**Aupcon Medical Technology (Wuxi) Co., Ltd.** is a global supplier of sports protection and rehabilitation products.

**including kinesiology tape, adhesive bandages, sports tape, leg sleeves, wrist braces, knee braces, orthopedic support products, etc.**

Aupcon provides breathable, comfortable and adjustable support for the knee, neck, leg, back, wrist, shoulder, waist, arm and other parts.

We are committed to improving our customers' comfort, support and ease of use. Our products are more user-friendly, easy to use and affordable, and are trusted by doctors, elite athletes and patients from all walks of life.

We combine innovation, efficiency and reliability to meet the needs of global sports brands, distributors and retailers. Contact an orthopedic brace manufacturer today for a competitive solution!

Our professional sales team can quickly respond to quotes, orders and technical guidance, ensure timely shipment, and assist with customs clearance when needed.

**13,000**

Square feet

**80+**

Modern machines

**100+**

Skilled employees



Strict quality control



OEM and ODM  
design services





# CONTENT

- 01      Cervical Collar
- 02      Shoulder Brace
- 03      Arm Brare
- 04      Lumbar Brace
- 05      Knee Orthosis
- 06      Ankle Brace



A woman with blonde hair tied back is sitting in a meditative lotus position on a grassy field. She is wearing a black sports bra with pink trim. Her eyes are closed, and her expression is peaceful. The background is a lush green park with tall trees and sunlight filtering through the leaves. A semi-transparent white banner with the text "CERVICAL COLLAR" in blue capital letters is positioned across the lower part of the image, partially covering her torso and the grass.

**CERVICAL COLLAR**



TJ-101

塑料颈托 II 型

Adjustable Cervical Collar



## Indication:

1. For the prevention and adjuvant treatment of cervical spondylosis, cervical muscle strain, such as cervical syndrome.
2. Cervical spine fracture fixation.
3. Cervical spondylosis surgery after rehabilitation period of the neck part fixing.

## Characteristic:

1. High-grade adjustment button can adjust the height of the neck, in accordance with the needs of different bone proportion of the population.
2. Back support board air holes design, improved greater lengthways vertical support.

**Size: Universal Size**

TJ-102

泡沫颈托 I 型

Philadelphia Cervical Collar



## Indication:

It is suitable for cervical fracture fixation, dislocation reduction and fixation, etc.

## Characteristic:

1. Imported British Foam (PLAZOTE).
2. 18 specifications meet the needs of different people.
3. Environmental material and Safe and it can direct touch the skin, nontoxic, no irritating smell, not containing CFC(chlorofluorocarbons) and HCFC(halogenated fluoride).
4. Without metal design, X-light, MRI and CT penetration effect are excellent.



TJ-102  
泡沫颈托 I 型  
Philadelphia Cervical Collar



Size:

型号(Model)	高度(Height)(cm)	颈围(Size ) (cm)
婴儿型(Infant)		
儿童型(Child)成人型(Adult)		
S	21/4"(5.7)	10"-13" (25-33)
S	31/4"(8.3)	10"-13"(25-33)
S	41/4"(10.8)	10"-13"(25-33)
S	51/4"(13.3)	10"-13"(25-33)
M	21/4"(5.7)	13" -16"(33-41)
M	31/4"(8.3)	13" -16"(33-41)
M	41/4"(10.8)	13"-16"(33-41)
M	51/4"(13.3)	13"-16"(33-41)
L	21/4"(5.7)	16"-18"(41-46)
L	31/4"(8.3)	16"-18"-(41-46)
L	41/4"(10.8)	16"-18" (41-46)
L	51/4"(13.3)	16"-18"(41-46)
XL	21/4"(5.7)	18"-19"(46-48)
XL	31/4"(8.3)	18"-19"(46-48)
XL	41/4"(10.8)	18"-19"(46-48)
XL	51/4"(13.3)	18"-19" (46-48)

TJ-103

塑料颈托 V 型

Cervical Collar



## Indication:

Ideal use for prevention and auxiliary treatment of cervical vertebra syndrome , such as cervical spondylosis and cervical muscle strain etc and neck support for cervical spondylosis rehabilitation after surgery.

## Characteristic:

1. Well ventilated support plate to increase cool and comfort.
2. Double velcro structure to ensure the best fixed effect.

Model	Size(cm)
S	37-42
M	42-47
L	47-52
X L	52-57



# TJ-103 (1)

## 塑料颈托 I 型

### Cervical Collar



### Indication:

1. For the prevention and adjuvant treatment of cervical spondylosis, cervical muscle strain, such as cervical syndrome.
2. Cervical spine fracture fixation.
3. Cervical spondylosis surgery after rehabilitation period of the neck fixing.

### Characteristic:

The ergonomic design of human body ensures the stability of the neck .

Soft sponge lining, good air permeability, comfortable to wear.

TJ-104

塑料颈托Ⅲ型

Emergency Cervical Collar



## Indication:

- Apply to cervical fracture fixation, dislocation reduction and fixation, etc

## Characteristic:

1. Four gears to adjust height
2. Without metal design, X- light, MRI and CT penetration effect are excellent.
3. Big trachea opening is convenient for carotid artery monitoring, and the open hole design in back is breathable and easy to diagnose.

**Size:** Universal size



## TJ-105 (2)

### 海绵颈托 I 型

### Soft Collar



### Indications:

- Ideal use for prevention and auxiliary treatment of cervical vertebra syndrome, such as cervical spondylosis and cervical muscle strain etc and neck support for cervical spondylosis rehabilitation after surgery.

### Characteristic:

- Designed according to the natural curve of head and neck, it is comfortable and easy to clean.

Size (规格)	Height(cm)	Length(cm)
Child S	5.5	22-24
Child M	6.5	24-26
Child L	7.5	26-28
Adult S	7.5	27-29
Adult M	8.5	29-33.7
Adult L	9	33.7-42
Adult XL	10.5	42-48

# TJ-105 (3)

## 海绵颈托 II 型

### Soft Collar



### Indications:

-Ideal use for prevention and auxiliary treatment of cervical vertebra syndrome, such as cervical spondylosis and cervical muscle strain etc and neck support for cervical spondylosis rehabilitation after surgery.

### Characteristic:

-Designed according to the natural curve of head and neck, it is comfortable and easy to clean.

Size (规格)	Height(cm)	Length(cm)
Child S	5.5	22-24
Child M	6.5	24-26
Child L	7.5	26-28
Adult S	7.5	27-29
Adult M	8.5	29-33.7
Adult L	9	37-42
Adult XL	10.5	42-48



TJ-106

颈托 II 型 (泡沫型)

Cervical Collar II (Inflatable)



## Indications:

- Suitable for cervical spondylosis, cervical vertebra bone hyperplasia and other people.

## Characteristic:

- According to ergonomic design the former is low, to avoid the extrusion of the neck sensitive parts, the use of safety, wearing more comfortable.
- Using large latex ball, inflatable time is short, streamlined design, beautiful and durable.

**Size:** Adult Universal size

TJ-107

## 塑料颈托Ⅳ型

Cervical Collar With Chin Support(Reinforced)



### Indications:

- Ideal use for prevention and auxiliary treatment of cervical vertebra syndrome , such as cervical spondylosis and cervical muscle strain etc and neck support for cervical spondylosis rehabilitation after surgery.

### Characteristic:

1. Well ventilated support plate to increase breathable and comfort.
2. Two pieces of self-adhesive support plate can adjust height freely to ensure perfect fit between chin and cervical collar.
3. Double velcro structure ensures the best fixation effect.

size	Neck Circumference (cm)
S	37-42
M	42-47
L	47-52
XL	52-57
XXL	57-62

**TJ-108**

**头颈胸固定器 I 型**

**Cervical Thoracic Orthosis I**



## **Indications:**

Use for patients with cervical spine injury and high thoracic injury, cervical spine (except for the atlantoaxial) injury after consolidation and rehabilitation of the campaign to protect, cervical arthritis and stability of the cervical spine fracture fixing and protecting treatment.

## **Characteristic:**

- 1.Product material: It made of the hardness of polymer composite materials.
- 2.The head ring fixation device to effectively limit the head rotation.
- 3.The chest compression height can be adjusted according to the patient height of the back side of the chest.
- 4.Occipital bone support height can adjust the height of occipital bone until it supports well.
- 5.Can be removed and thickened breathable lining, make wearing more comfortable, and can be cleaned.

**Size:** Adult M、L



TJ-109

头颈胸固定器 II 型

Cervical Thoracic Orthosis II



## Indications:

Use for patients with cervical spine injury and high thoracic injury, cervical spine (except for the atlantoaxial) injury after consolidation and rehabilitation of the campaign to protect, cervical arthritis and stability of the cervical spine fracture fixing and protecting treatment.

## Characteristic:

- 1.Product material: It made of the hardness of polymer composite materials.
- 2.The head ring fixation device to effectively limit the head rotation.
- 3.The chest compression height can be adjusted according to the patient height of the back side of the chest.
- 4.Occipital bone support height can adjust the height of occipital bone until it supports well.
- 5.Can be removed and thickened breathable lining, make wearing more comfortable, and can be cleaned.

**Size:** Adult、Child



**SHOULDER BRACE**



TJ-201

臂外展固定器 I - I 型

Humeral Abduction Orthosis



## Indications:

- Deltoid muscle paralysis, supraspinous tendon rupture, shoulder joint fracture and dislocation reduction, shoulder joint surgery, brachial plexus injury patients.

## Characteristic:

1. Abduction angle and elbow flexion angle adjustable through the universal wheel.
2. Liner cotton pad can remove and wash.
3. Suitable for left and right hands.

**Size:** Adult Universal size



TJ-201 (1)

臂外展固定器 I - II 型

Humeral Abduction Orthosis



## Indications:

- It is suitable for the maintenance of functional position in patients with deltoid palsy, supraspinal tendon fracture, shoulder joint fracture, dislocation repair, shoulder joint surgery and brachial plexus nerve injury.

## Characteristic:

- Soft composite material, comfortable to wear, breathable and skinfriendly.
- Velcro is strong and durable with long service life.
- The soft pad for hands is made of skin-friendly breathable fabric to prevent hand soreness after long time placement.
- Triangle stable support, angle can be exchanged

**Size:** Adult Universal size

TJ-202

臂外展固定器 II - I 型

Shoulder Abduction Orthosis



## Indications:

- It is suitable for the fixation of the shoulder and the shoulder joint dislocation after the reduction and fixation of the shoulder joint and the upper end of the humerus.

## Characteristic:

- Composite sponge pillow fixed shoulderjoint outside booth,so general, sticky connection can easily adjust the abduction angle. Use the ball, promote limb activity,prevent muscle atrophy.

**Size:** Universal size

TJ-202 (2)

臂外展固定器 II - II 型

Shoulder Abduction Orthosis



### **Indications:**

- It is suitable for the fixation of the shoulder and the shoulder joint dislocation after the reduction and fixation of the shoulder joint and the upper end of the humerus.

### **Characteristic:**

- Composite sponge pillow fixed shoulder joint outside booth, so general, sticky connection can easily adjust the abduction angle. Use the ball, promote limb activity, prevent muscle atrophy.

**Size:** S、M、L、XL



TJ-202 (3)

臂外展固定器 II-III 型

Shoulder Abduction Orthosis



## Indications:

- Suitable for fixation after rotator cuff repair, dislocation of shoulder joint, reduction and fracture of shoulder and upper humerus.

## Characteristic:

- Composite sponge pillow for shoulder joint external placement.
- The design is in line with the physiological characteristics of the human body and is comfortable and comfortable. The sleeper is fastened by the shoulder strap to avoid slipping.
- It is comfortable and portable and easy to use.

**Size:** Left/Right

TJ-203

肩部固定带 I 型

Shoulder Brace



## Indications:

- Suitable for the shoulder pain and chronic pain, inflammation, shoulder joint injury, shoulder joint instability or habitual shoulder joint .

## Characteristic:

- The shoulder and the upper arm of the two adjustable band, can adjust the size of the right, left and right side of the design, so that more convenient to wear

**Size:** Universal size

# TJ-203 (1)

## 肩部固定带 II 型

### Shoulder Brace



## Indications:

- Postural control in patients with mild hunchback.
- Control of correct posture in developing adolescents.

## Characteristic:

- Light and breathable, integral fit, can be worn in clothes, strong concealment, easy to wear.

(Size)	Chest circumference(CM)
XS	55-65
S	65-80
M	80-95
L	95-115



TJ-204

锁骨固定带Ⅲ型

Clavicle Brace



## Indications:

- Clavicle and scapula fracture, arched back and back or hunchback.

## Characteristic:

- Use high quality composite flannelette upset paragraph, various molding, metal ring, which divides the straps, convenient adjustment, to assist in the fixation of fracture of clavicle also remind posture.

**Size:** S、M、L、XL

TJ-204 (1)  
锁骨固定带IV型  
Clavicle Brace



**Indications:**

- Clavicle and scapula fracture, arched back and back or hunchback.

**Characteristic:**

- Use high quality composite flannelette upset paragraph, various molding, metal ring, which divides the straps, convenient adjustment, to assist in the fixation of fracture of clavicle also remind posture.

**Size:** Universal size

TJ-205 (2)

背部固定器Ⅲ型

Spine Orthopedics Brace



## Indications:

- Postural control in patients with mild humpback.

Correct posture control in adolescence.

## Characteristic:

- lightweight and breathable, with overall fit. It can be worn in clothes, concealed and easy to wear.

**Size:** S、M、L、XL、XXL





**ARM BRACE**

TJ-301

可调式肘关节固定器 I 型

Adjustable Elbow Immobilizer



## Indication:

1. Conservative treatment of unstable elbow joint with or without other injuries.
2. Post-operative care for elbow joint.
3. Functional follow-up treatment of stable fractures.
4. Post-operative care for ligament instability of elbow joint.

## Characteristic:

- Individual early functional treatment of the activity and correction.
- Fixed elbow joint by fixing the upper arm and forearm. The elbow is limited to the range of 0 to 120 degrees. Ensure the flexibility of in the early stage. Enhance elbow joint activity to speed up the recovery.
- Adjust the joint by simply pressing the button.



Size \_\_\_\_\_

Universal size \_\_\_\_\_



TJ-301 (1)

可调式肘关节固定器 II 型

Adjustable Elbow Immobilizer



## Indication:

- 1.Conservative treatment of elbow instability with or without other injuries.
- 2.Post-operative care of elbow joint stabilization.
- 3.Functional follow-up treatment of stable fractures.
- 4.Post-operative care of elbow ligament instability.

## Characteristic:

- 1.Individual early functional therapy activities and orthopedics
- 2.Fixed elbow joint by fixing the upper arm and forearm.The elbow is limited to the range of 0 to 120 degrees.Ensure the flexibility of in the early stage. Enhance elbow joint activity to speed up the recovery.

Size \_\_\_\_\_

Universal size ( Left/Right)



TJ-301 (2)

可调式肘关节固定器Ⅲ型

Adjustable Elbow Immobilizer



## Indication:

1. Conservative treatment of unstable elbow joint with or without other injuries.
2. Post-operative care for elbow joint.
3. Functional follow-up treatment of stable fractures.
4. Post-operative care for ligament instability of elbow joint.

## Characteristic:

- Individual early functional treatment of the activity and correction.
- Fixed elbow joint by fixing the upper arm and forearm. The elbow is limited to the range of 0 to 120 degrees. Ensure the flexibility of in the early stage. Enhance elbow joint activity to speed up the recovery.

Size \_\_\_\_\_

Universal size ( Left/Right )

## TJ-302 (3) 腕关节固定带 I 型 Wrist Brace



### Indications:

- Post-traumatic or post-operative wrist Chronic soft tissue injury;
- Wrist rheumatoid arthritis;
- Solid after removing gypsum bandage.

### Characteristic:

- Support and fix wrist joint;
- Three extra-wide support belts placed on the wrist strengthen the fixation effect;
- Plastic aluminum plate ensures proper fixed position.

Size	Wrist circumference(cm)
S (Left/Right)	13-15
M(Left/Right)	15-18
L (Left/Right)	18-21
XL (Left/Right)	>21

TJ-302 (4)

腕关节固定带 II 型

Wrist Brace



## Indications:

- Chronic wrist after or after trauma;
- Soft tissue injury;
- Wrist rheumatoid arthritis;
- Fixation after removing the plaster bandage.

## Characteristic:

- It is easy to wear, and the rope pulling device improves wrist stability and patient comfort.

<u>Size</u>	<u>Wrist circumference(cm)</u>
<u>Child Universal Size</u>	<u>12-15</u>
<u>Adult Universal Size</u>	<u>16-23</u>



# TJ-302 (S)

## 腕关节固定带VI型

### Wrist Brace



### Indications:

- Chronic soft tissue injury of wrist after or after trauma;
- Wrist rheumatoid arthritis;
- Fixation after removing gypsum bandage

### Characteristic:

- Skin-friendly fabric, soft and comfortable;
- Removable aluminum plate ensures proper fixed position;
- The drawstring design is convenient to adjust.

<u>Size</u>	<u>Wrist circumference(cm)</u>
S (Left/Right)	<13
M(Left/Right)	13-25
L (Left/Right)	>24.5

TJ-303

腕关节固定带Ⅷ型

Adjustable Forearm Support (Enhanced)



## Indications:

·Nursing for wrist fracture, dislocation of metacarpal bone and ligament injury and rehabilitaion.

## Characteristic:

·Forearm protection part of the extension to ensure the rehabilitation process.

·Aluminum alloy arc support in accordance with the hands of the curve, wearing more comfortable.

Size	Wristcircumference(cm)
S (Left/Right)	13-15
M(Left/Right)	15-18
L (Left/Right)	18-21
XL (Left/Right)	>21

TJ-304 (1)

拇指固定带 II 型

Thumb Protect Brace



## Indications:

- Sprained thumb joint fixed
- Tenosynovitis prevention and protection
- Fixation of sports injury
- Mouse hand prevention and protection

## Characteristic:

- Breathable fabric, soft and skin-friendly.
- Injection molded board with ergonomic design, more comfortable to wear.
- Cyclotron wrist strap for better wrist fixation.
- Humanized thumb strap, easy to fix thumb, avoid tendon sheath friction.

Size \_\_\_\_\_

Universal size \_\_\_\_\_



## TJ-304 (2) 拇指固定带Ⅲ型 Thumb Brace



### **Indications:**

-Pain in thumb root, mild thumb sprain, tenosynovitis, postoperative fixation.

### **Characteristic:**

- The opening design at the thumb can prevent pressure on the painful part;
- After wearing it, it does not affect the holding action, ensuring the range of wrist movement;
- It is light and breathable, and can effectively sweat.

Size

Universal size(Left/Right)

TJ-305 (1)

腕关节固定带IV型

Thumb Protect Brace(Reinforced)



## Indications:

-Thumb sprain, arthritis, tenosynovitis, fasciitis, stable fracture and joint release.

## Characteristic:

- Knob type fine-tuning tightening system, easy to wear, uniform pressure, easy to adjust the tightness.
- The aluminum support conforms to the physiological design, and it is tightly fitted and fixed.
- High quality fabric, skin friendly lining, moisture absorption and sweat, comfortable and breathable
- Fine workmanship, light and short, stylish appearance.

Size \_\_\_\_\_

S (Left/Right) \_\_\_\_\_

L(Left/Right) \_\_\_\_\_

TJ-305 (2)

腕关节固定带 V 型

Thumb Brace



## Indications:

- Chronic soft tissue injury of wrist after trauma or surgery;
- Rheumatoid arthritis of wrist;
- Fixation after removing plaster bandage.

## Characteristic:

- Skin-friendly fabric, soft and comfortable
- Drawstring design is convenient to adjust.

Size Wristcircumference(cm)

S (Left/Right) <13

M(Left/Right) 13-25

L (Left/Right) >24.5

TJ-306

前臂固定带VI型

Multi-function Arm sling (Sponge Type)



## Indications:

- Shoulder contusion and strain.
- Shoulder dislocation.
- The stability of the top of the humerus fracture. Fracture of the shoulder blade.

## Characteristic:

- Firm fixation of shoulder or upper arm.
- Chest fixed strap to strengthen the stability. Make sure that no affect to checking wound.

Size \_\_\_\_\_

S、 M、 L、 XL



TJ-306 (1)

前臂固定带VII型

Multi-function Arm sling (Sponge Type)



## Indications:

- Shoulder contusion and strain.
- Shoulder dislocation.
- The stability of the top of the humerus fracture.Fracture of the shoulder blade.

## Characteristic:

- Firm fixation of shoulder or upper arm.
- Chest fixed strap to strengthen the stability. Make sure that no affect to checking wound.

Size \_\_\_\_\_

Universal size \_\_\_\_\_

TJ-307 (2)

前臂固定带Ⅷ型

Arm Sling (Mesh)



## Indications:

-Shoulder dislocation fixation, arm fracture fixation, wrist sprain fixation, sports sprain fixation。

## Characteristic:

-Skin-friendly fabric, light and breathable, separated shoulder sleeve, evenly distribute weight, reduce pain and fatigue, the length of front and back shoulder straps can be adjusted, comfortable and fixed, fit both left and right.

Size \_\_\_\_\_

Universal size \_\_\_\_\_

TJ-307 (3)

前臂固定带Ⅲ型

Arm Sling (Mesh)



## Indications:

- Rotator cuff reconstruction.
- Fixation of shoulder joint dislocation.

## Characteristic:

- The main force goes to the shoulder.
- Avoid neck fatigue caused by concentrated pressure on neck.
- It can be adjusted according to the angle of elbow joint and the position of the hand.

Size \_\_\_\_\_

Universal size \_\_\_\_\_

## TJ-308 (1) 前臂固定带IV型 Arm Sling



### Indications:

-Shoulder dislocation fixation, arm fracture fixation, wrist sprain fixation, sports sprain fixation, etc.

### Characteristic:

-Skin-friendly fabric, light and breathable, adjustable length of front and back shoulder straps, comfortable and fixed, available for both hand.

Size

Universal size



## TJ-308 (2) 前臂固定带X型 Arm Sling



### **Indications:**

-Shoulder dislocation fixation, arm fracture fixation, wrist sprain fixation, sports sprain fixation, etc.

### **Characteristic:**

-Skin-friendly fabric, light and breathable, comfortable and fixed, available for both hand.

Size \_\_\_\_\_

Child Universal Size

TJ-309

前臂固定带 V 型

Triangle Arm Sling



## Indications:

- The forearm fractures
- Wrist fracture
- Wrist sprain and finger sprain

## Characteristic:

- Adjust according to elbow Angle
- Can replace triangle bandage, easy to use, simple and comfortable.

Size \_\_\_\_\_

Child Universal Size \_\_\_\_\_

Adult Universal Size \_\_\_\_\_





**LUMBAR BRACE**



TJ-401 (1)

腰部固定器 I 型

Lumbar Support Orthosis( I )



## Indications:

-Lumbar spine protection, lumbar postoperative support fixation, disc degeneration disease, disc swelling, fracture treatment, osteoporosis, spinal injury, kyphosis, spinal loosening, spondylolisthesis, vertebral forward movement, spinal stenosis, vertebral disease (osteoarthritis/degenerative joint disease)

## Characteristic:

- This lumbosacral integral fixator is designed for control, stability, and severe pain relief related to common spinal problems.
- This product adopts pulley rope design, effectively reduce friction, more flexible and stable operation.
- The rear panel can be raised or lowered through the knob, the height is adjustable, convenient and fast.
- The inner liner adopts hot pressing technology, which is ergonomic design, making it more comfortable to wear.
- Code design, adapt to a wide crowd.

Size	Waist circumference (cm)
Universal Size	70-120



TJ-401 (2)

拉绳盆骨固定器 I 型

Lumbar Support Orthosis ( I )



## Indications:

-Lumbago sour, pubic separation, pubic pain, hip width, etc.

## Characteristic:

·Use high quality material for pelvic support, easy to wear.

·Pulley rope design, effectively reduce the friction, more flexible and stable operation, wearing more comfortable.



Size	Hip circumference(cm)
M	80-110
L	110-130

TJ-402

腰部固定器 II 型

Lumbar Support Orthosis ( II )



## Indications:

-Lumbar protection, lumbar post-operative support and fixation, degenerative disc disease, prolapse of intervertebral disc, fracture management, osteoporosis, spinal injury, kyphosis, spondylolysis, spondylolisthesis, spondylolisthesis, spinal stenosis, spondylopathy (osteoarthritis/degenerative joint disease)

## Characteristic:

- This lumbosacral integral fixator is designed for control, stability, and severe pain relief related to common spinal problems.
- This product adopts pulley rope design to control the whole waist and back support plate can be disassembled.

Model	Waist circumference (cm)
XS	50-65
S	65-80
M	80-95
L	95-110
XL	110-125
XXL	>125

TJ-402 (2)

腰背固定器 II 型

Lumbar Support Orthosis( II )



## Indications:

- Juvenile compulsory spondylitis As(rickets)
- Mechanical pain in back
- PVP minimally invasive surgery
- Osteoporosis

## Characteristic:

- The rope pulling device improves waist stability and patient comfort.
- Optimize PVP effect



Model	Waist circumference (cm)
XS	50-65
S	65-80
M	80-95
L	95-110
XL	110-125
XXL	>125



TJ-403

腰部固定器IV型

Lumbar Support Orthosis (IV)



## Indications:

- Early osteoporosis, compression fracture of thoracic vertebra kyphosis and thoracic mechanical pain.
- Adjunctive treatment for pre-operation and post-operation.

## Characteristic:

- This product can relieve pain by hyperextension of chest.
- The rear plastic plate can effectively release the spine and can correct 15 degrees of spinal lordosis.
- The correct range is from shoulder to the lower spine and then sacrum.



Model	Waist circumference (cm)
XS	50-65
S	65-80
M	80-95
L	95-110
XL	110-125
XXL	>125



TJ-403 (1)

腰部固定器VII型

Lumbar Support Orthosis (VII)



## Indications:

-Lumbar protection, lumbar post-operative support and fixation, degenerative disc disease, prolapse of intervertebral disc, fracture management, osteoporosis, spinal injury, kyphosis, spondylolysis, spondylolisthesis, spondylolisthesis, spinal stenosis, spondylopathy (osteoarthritis/degenerative joint disease).

## Characteristic:

- Comfortable to wear, the effect is prominent according to the body needs to adjust the belt binding force size, more comfortable than the steel belt.
- Ergonomic design, according to human lumbar convex design, closely fit lumbar lordosis, perfect fit.
- Pull rope design pressure injection molding support plate forward, provide super support.
- Smooth mechanical pulley contraction system, convenient for users to flexibly adjust the tension through the rope

Model	Waist circumference (cm)
XS	50-65
S	65-80
M	80-95
L	95-110
XL	110-125
XXL	>125

TJ-404

胸腰椎固定器 I 型

Thoracolumbar Sacral Orthosis( I )



## Indications:

-Thoracolumbar post-operation, compression fracture,degenerative lumbar disc herniation, congenital disease and other thoracic lumbar problems.

## Characteristic:

-The humanity design remedies the lack of customizedlumbar support by human body curve.Provide morecomfortable and convenient service for all patients who has lumbar problems.

-It is designated as male and female type.

Size	Waist	Chest	Front Height	Back Height	Unit
S	45-60	60	24	26	cm
M	55-70	70	28	30	cm
L	65-80	80	33	35	cm

Size	Waist	Chest	Front Height	Back Height	Unit
S	71-78	74	40	46	cm
M	78-86	82	42	48	cm
L	86-94	90	42	48	cm
XL	94-102	100	46	50	cm

TJ-404 (1)

胸腰椎固定器Ⅲ型

Thoracolumbar Sacral Orthosis(Ⅲ)



## Indications:

-Thoracolumbar post-operation, compression fracture, degenerative lumbar disc herniation, congenital disease and other thoracic lumbar problems.

## Characteristic:

- Built-in integral air bag, manual air charging and discharging for the lumbar support more fitable, skin friendly breathable lining, make wearing more comfortable.
- Two sides of three groups of velcro fixed, convenient and firm.



Size \_\_\_\_\_

Universal size \_\_\_\_\_



TJ-404 (2)

胸腰椎固定器 V 型

Thoracolumbar Sacral Orthosis(V)



## Indications

- Thoracolumbar post-operation, compression fracture, degenerative lumbar disc herniation, congenital disease and other thoracic lumbar problems.

## Characteristic

- Built-in integral air bag, manual air charging and discharging, for the lumbar support more fitable, skin friendly breathable lining, make wearing more comfortable.
- Two sides of three groups of velcro fixed, convenient and firm.



## Size

Universal Size



TJ-404 (3)

胸腰椎固定器 VI 型

Thoracolumbar Sacral Orthosis(VI)



## Indications

- Thoracolumbar post-operation, compression fracture, degenerative lumbar disc herniation, congenital disease and other thoracic lumbar problems.

## Characteristic

- Humanized design provides more convenient service for all thoracic and lumbar spine patients.
- Ergonomic design, front adjusting hinge makes wearing more comfortable.



## Size

Universal Size

TJ-405

胸腰椎固定器 II 型

Lumbar Support Orthosis( II )



## Indications

- Thoracolumbar post-operation, compression fracture, degenerative lumbar disc herniation, congenital disease and other thoracic lumbar problems.

## Characteristic

- Humanized design provides more convenient service for all thoracic and lumbar spine patients.

## Size

Size	Waist	Chest	Front Height	Back Height	Unit
S	66-74	74-81	25	30	cm
M	74-81	81-89	25	30	cm
L	81-89	89-96	25	30	cm
XL	89-96	96-104	25	30	cm



TJ-405 (1)

胸腰椎固定器VII型

Lumbar Support Orthosis(VII)



## Indications

- Thoracolumbar post-operation, compression fracture, degenerative lumbar disc herniation, congenital disease and other thoracic lumbar problems.

## Characteristic

- Built-in integral air bag, manual air charging and discharging, for the lumbar support more fitable higher, skin friendly breathable lining, make wearing more comfortable.
- Two sides of three groups of Velcro fixed, convenient and firm.



## Size

Universal Size

TJ-406

胸腰椎固定器Ⅷ型

Lumbar Support



## Indications

- For thoracolumbar vertebral fracture, osteoporosis and vertebral displacement.

## Characteristic

- The formed hard plate is made in accordance with ergonomic design, which provides a stable support for the lumbar sacral and abdomen.



## Size

Size	Waist circumference (cm)
M	60-90
L	85-115



TJ-406 (1)

胸腰椎固定器IX型

Lumbar Support



## Indications

- Lumbar spondylolisthesis, lumbar disc herniation, low back pain, muscle dysfunction, lumbar joint pain and sports protection.

## Characteristic

- Wear comfortable, strong support, can be adjusted by both sides of the velcro and adjusting rod comfort, high-end product design.



## Size

Size	Waist circumference(cm)
M	60-90
L	85-115

TJ-407

腹带XIII型

Ostomy Abdominal Brace



## Indications

- It is for colostomy patients to wear on the abdominal stomato fix bag and prevent parastomal hernia and hernia.

## Characteristic:

- It is flexible and can be used for different body size.
- Unisex.

## Size

Size	Waist circumference(cm)
XS	<85
S	85-95
M	95-105
L	105-115
XL	115-125

**TJ-408**  
**腹带V型**  
**Abdominal Brace**



## Indications:

- Help support abdominal weight, lose weight and trim belly.

## Characteristic:

- Support weight of the abdomen and provide stability of the back. The bilateral scalable sticky buckle can provide additional length.

## Size

Size	Waist circumference(cm)
S	65-75
M	75-85
L	85-95
XL	95-105
XXL	105-115

TJ-408 (3)  
腹帶XI型  
Abdominal Brace



**Indications:**

- Postpartum and postoperative abdominal protrusion, compression of the abdomen.

**Characteristic:**

- Wide design with detachable auxiliary pad. It is natural and comfortable. It makes patients feel breathable. It is elastic.
- Elastic support strips provide proper support. They are not easy to deform. High-quality fish ribbon fabric improves the breathability of the product.

**Size**

Size	Abdominal Girth (cm)
S	60-90
M	90-120
L	120-150



TJ-408 (4)

腹带XII型

Abdominal Brace



## Indications:

- Protection of surgical or postpartum wounds, compression of the abdomen.

## Characteristic:

- It supports the weight of the abdomen and provides back stability, It is made by high elastic fish ribbon, a wide range of adjustment, more comfortable to wear.

## Size

Size	Abdominal Girth (cm)
S	60-90
M	90-120
L	120-150

TJ-409

腰椎固定带 I - II 型

High Back Brace



## Indications:

- Lumbar soft tissue injury
- Lumbar disc herniation
- Lower thoracic spinal injury or slippage
- Mild skeletal injury rehabilitation.

## Characteristic:

- The higher design can enlarge support range to lower thoracic vertebrae.

## Size

Size	Waist circumference(cm)
S	75-85
M	85-95
L	95-105
XL	105-115
XXL	115-125

TJ-411

腰椎固定带 I - I 型

Back Brace



## Indications:

- Injuries of lumbar and sacral vertebrae, such as acute or chronic low back pain, preoperative or postoperative treatment. The pressure can be adjusted.

## Characteristic:

- This product is suitable for patients with strong body.

## Size

Size	Waist circumference(cm)
S	75-85
M	85-95
L	95-105
XL	105-115
XXL	115-125
XXXL	125-135

TJ-413

腰椎固定带 I-VI型

Back Support Belt(Elastic)



## Indications:

- Low back pain
- muscle dysfunction
- Lumbar spondylolisthesis
- lumbar intervertebral disc disease(herniation, intervertebral disc herniation)
- Lumbar and sacral joint pain.

## Characteristic:

- The waistline is ergonomic for better support.
- There are sponges on the back to cushion back pressure.

## Size

Size	Waist circumference(cm)
S	75-85
M	85-95
L	95-105
XL	105-115
XXL	115-125
XXXL	125-135



TJ-413 (1)

腰椎固定带 I-IV 型

Back Support Belt(Elastic)



## Indications:

- Low back pain, muscle dysfunction
- lumbar spondylolisthesis,
- lumbar intervertebral disc disease(herniation, intervertebral disc herniation)
- lumbar and sacral joint pain and sports protection.

## Characteristic:

- Integrated cutting, 4 aluminum stays + 2 PP stays, evenly share the waist pressure, cross type reinforcement belt, according to the need to adjust the pressure, fixed support, using high-quality breathable material, wearing is not stuffy.

## Size

Size	Waist circumference(cm)
S	75-85
M	85-95
L	95-105
XL	105-115

TJ-414

腰椎固定带 I-V 型

Back Brace



## Indications:

- Low back pain, muscle dysfunction, lumbar spondylolisthesis, lumbar intervertebral disc disease (herniation, intervertebral disc herniation), lumbar and sacral joint pain.

## Characteristic:

- Strong support strips
- Breathable pocket is used to do hot or cold compress therapy.

## Size

Size	Waist circumference(cm)
S	75-85
M	85-95
L	95-105
XL	105-115
XXL	115-125
XXXL	125-135
XXXXL	>135

TJ-416

腰椎固定带 I—IX型

Waist Support Belt (Reinforced)



## Indications:

- Lumbar pain, long-term excessive force protection.

## Characteristic:

- This product is used in the work of safety protection, with strong flexibility, more suitable for manual workers.

## Size

Size	Waist circumference(cm)
S	70-80
M	80-90
L	90-100
XL	100-110
XXL	110-120

TJ-417 (2)

肋骨固定帶Ⅲ型

Rib BraceⅢ



## Indication:

- Fixed rib injury.
- Damage to chest tissue.
- Postoperative fixation of fractures.

## Characteristic:

- Convenient to use, comfortable to wear, large adjustment range.

## Size

S、M、L、XL



TJ-418 (1)

腰椎固定带 II—II 型

Lumbar Brace



## Indication:

- Low back pain, muscle insufficiency, lumbar spondylolisthesis, intervertebral disc disease, fractures, lumbar joint pain, etc.

## Characteristic:

- The heat-sealing process does not allow balls, slips, and is firmer.

The high-stretch woven fabric has stronger stretching force and is comfortable to wear.

## Size

Size	Waist circumference(cm)
S	57-67
M	80-90
L	90-100
XL	100-110
XXL	110-120



## KNEE ORTHOSIS



TJ-501

可调式膝关节固定器IV型

Adjustable Knee Orthosis



## Indication:

- Conservative treatment of comprehensive knee instability.
- Repair and fixation of fractures around the knee joint, soft tissue injuries, skeletal fractures, dislocations, meniscus injuries, muscle contractures, arthritis, knee flexion and other diseases.

## Characteristic:

- Help the knee joint to restore the original function or functional exercise during the rehabilitation period.
- The length can be adjusted according to personal needs.
- The adjustable chucks on both sides can move freely from  $0^{\circ}$  - $120^{\circ}$  , and the angle can be adjusted and fixed within the range of  $0^{\circ}$  - $90^{\circ}$  .

## Size

Universal size



TJ-S01 (2)

可调式膝关节固定器VI型

Adjustable Knee Orthosis



## Indication:

- Conservative treatment of comprehensive knee instability.
- Repair and fixation of fractures around the knee joint, soft tissue injuries, skeletal fractures, dislocations, meniscus injuries, muscle contractures, arthritis, knee flexion and other diseases.

## Characteristic:

- Help the knee joint to restore the original function or functional exercise during the rehabilitation period.
- The length can be adjusted according to personal needs.
- The adjustable chucks on both sides can move freely from  $0^{\circ}$  - $120^{\circ}$  , and the angle can be adjusted and fixed within the range of  $0^{\circ}$  - $90^{\circ}$  .

## Size

Universal size





TJ-S01 (3)

可调式膝关节固定器Ⅶ型

Adjustable Knee Orthosis



## Indication:

- Conservative treatment of comprehensive knee instability.
- Repair and fixation of fractures around the knee joint, soft tissue injuries, skeletal fractures, dislocations, meniscus injuries, muscle contractures, arthritis, knee flexion and other diseases.

## Characteristic:

- Help the knee joint to restore the original function or functional exercise during the rehabilitation period.
- The length can be adjusted according to personal needs.
- The adjustable chucks on both sides can move freely from  $0^{\circ}$  - $120^{\circ}$  , and the angle can be adjusted and fixed within the range of  $0^{\circ}$  - $90^{\circ}$  .

## Size

Universal size



TJ-502

可调式膝关节固定器Ⅲ型

Adjustable Knee Orthosis



## Indication:

- Conservative treatment of comprehensive knee instability.
- Repair and fixation of fractures around the knee joint, soft tissue injuries, skeletal fractures, dislocations, meniscus injuries, muscle contractures, arthritis, knee flexion and other diseases.

## Characteristic:

- Help the knee joint to restore the original function or functional exercise during the rehabilitation period.
- The length can be adjusted according to personal needs.
- The adjustable chucks on both sides can move freely from  $0^{\circ}$  - $120^{\circ}$  , and the angle can be adjusted and fixed within the range of  $0^{\circ}$  - $90^{\circ}$  .

## Size

Universal Size



TJ-S02 (2)

可调式膝关节固定器IX型

Adjustable Knee Orthosis



## Indication:

- Repair and fixation of fractures around the knee joint, soft tissue injuries, skeletal fractures, dislocations, meniscus injuries, muscle contractures, arthritis, knee flexion and other diseases.

## Characteristic:

- Avoid hyperextension and limit flexion and extension.
- The long brackets on both sides can prevent the knee joint from turning in and out.
- The four-point principle can stabilize the knee joint due to the instability of the knee joint caused by osteoarthritis.

## Size

Universal Size



TJ-502 (3)

可调式膝关节固定器X型

Adjustable Knee Orthosis



## Indication:

- Conservative treatment of comprehensive knee instability.
- Repair and fixation of fractures around the knee joint, soft tissue injuries, skeletal fractures, dislocations, meniscus injuries, muscle contractures, arthritis, knee flexion and other diseases.

## Characteristic:

- Help the knee joint to restore the original function or functional exercise during the rehabilitation period.
- The length can be adjusted according to personal needs.
- The adjustable chucks on both sides can move freely from  $0^{\circ}$  - $120^{\circ}$  , and the angle can be adjusted and fixed within the range of  $0^{\circ}$  - $90^{\circ}$  .

## Size

Universal Size





TJ-S02 (4)

可调式膝关节固定器XIII型

Adjustable Knee Orthosis



## Indication:

- Conservative treatment of comprehensive knee instability.
- Repair and fixation of fractures around the knee joint, soft tissue injuries, skeletal fractures, dislocations, meniscus injuries, muscle contractures, arthritis, knee flexion and other diseases.

## Characteristic:

- Help the knee joint to restore the original function or functional exercise during the rehabilitation period.
- The length can be adjusted according to personal needs.
- Using carbon fiber support, high strength and low density, more light and stable
- The adjustable chucks on both sides can move freely from  $0^{\circ}$  - $120^{\circ}$  , and the angle can be adjusted and fixed within the range of  $0^{\circ}$  - $90^{\circ}$  .

## Size

Universal Size



TJ-S03

可调式膝关节固定器 I 型

Adjustable Knee Orthosis



## Indication:

- Repair and fixation of fractures around the knee joint, soft tissue injuries, skeletal fractures, dislocations, meniscus injuries, muscle contractures, arthritis, knee flexion and other diseases.

## Characteristic:

- Avoid hyperextension and limit flexion and extension.
- The long brackets on both sides can prevent the knee joint from turning in and out.
- The four-point principle can stabilize the knee joint due to the instability of the knee joint caused by osteoarthritis.

## Size:

M	41-47 cm
L	47-53 cm

TJ-505

可调式膝关节固定器XI型

Short Knee Orthosis



### Indication:

- Conservative treatment of collateral ligament injury.
- Meniscus injury.
- Traumatic or degenerative lesions of the knee joint.

### Characteristic:

- The knee joint is fixed on the side stably.
- Prevent excessive stretching.

### Size:

S、M、L、XL

TJ-507

膝部固定带 V 型

Kneepad



### **Indication:**

- Patella malacia.
- Joint edema, chronic inflammation after injury.
- The patella is moved up, subluxed or dislocated.

### **Characteristic:**

- Knee socket open design.
- There are live buckles and pressure tape on both sides of the middle to tighten the patella and knee.

### **Size:**

Universal size



TJ-508

膝部固定带 I 型

Knee Immobilizer



## Indication:

- Dislocation of the patella.
- Conservative treatment after ligament injury.
- Acute anterior knee pain.

## Characteristic:

- Reduce the tension of the posterior cruciate ligament (PCL) and reduce the pressure on the patella joint.
- The fixed branches on the inside, outside and back can be disassembled according to the disease.

## Size:

S 45cm

M 55cm

L 65cm

XL 75cm

TJ-508 (2)

膝部固定帶Ⅲ型

Knee Immobilizer



### Indication:

- Dislocation of the patella.
- Conservative treatment after ligament injury.
- Acute anterior knee pain.

### Characteristic:

- Reduce the tension of the posterior cruciate ligament (PCL) and reduce the pressure on the patella joint.
- The fixed branches on the inside, outside and back can be disassembled according to the disease.

### Size:

S	40cm
M	50cm
L	60cm
XL	65cm

TJ-509

膝踝足固定器

Knee Ankle-foot Orthosis



## Indication:

- Acute immobilization and protection of knee and tibiofibular injuries;
- Conservative treatment of the knee and tibiofibular injuries, soft tissue injuries, nerve injuries, crush injuries, or functional external fixation after surgery.

## Characteristic:

- Made of medium-hardness polymer material, the fixing effect is reliable;
- Thick removable and washable lining makes wearing more comfortable;
- Adjustable angle;
- Wearing it will not affect X-ray,CT and MRI examinations.

## Size:

Light/Right

# TJ-S10

## 髌外展固定器（单髌）

## Single Hip Orthosis



### Indications:

- Hip injury, soft tissue injury of the hip joint;
- Protective fixation after hip replacement, prevent displacement or dislocation of hip fractures.

### Characteristic:

- Using positioning disc locking type hip ammonium chain, which can fix the hip joint at multiple different flexion angles;
- The stainless steel support bar is equipped with a chute for easy length adjustment;
- Thickened lining makes it comfortable to wear and easy to disassemble;
- The waist and leg straps can be adjusted according to the patient's body shape

### Size:

Universal size



**TJ-S10 (1)**

**髋外展固定器（双髋）**

**Double Hip Orthosis**



### **Size:**

**Child:**S、 M、 L

**Adult universal size**

### **Indications:**

- Conservative fixation of hip joint and surrounding soft tissue injuries;
- Fixation after hip replacement;
- Early standing during rehabilitation training for patients with hip dislocation.

### **Characteristic:**

- Using positioning disc locking type hip ammonium chain, which can fix the hip joint at multiple different flexion angles;
- The stainless steel support bar is equipped with a chute for easy length adjustment;
- Thickened lining makes it comfortable to wear and easy to disassemble;
- The waist and leg straps can be adjusted according to the patient's body shape.





**ANKLE BRACE**



TJ-601 (1)

跟腱靴固定器VI型

Short Knee Orthosis



## Indication:

- Use for fixation of foot and ankle fracture, lower tibia and fibula fracture, ankle ligament damage, trauma of Achilles tendon fracture, stable fracture and postoperative rehabilitation etc.

## Characteristic:

- Adjustable button is easy to operate which can be locked in a certain angle or limited range and can be adjusted protected range-of-motion according to patient's condition of rehabilitation;
- Inflatable soft cushion is made of composite material, soft and comfortable. Filling and discharging operation is simple with the hidden inflatable air bag, which can be adjusted according to the needs of the patient's recovery. The fixed effect is very obvious
- Anti-slip design of sole, enhance the stability
- Universal design fits right and left foot

Size	Maximum foot length(cm)	Shoe size
S	24.5	36-40
M	26.5	40-42
L	29	42-45

TJ-602 (1)

跟腱靴固定器Ⅲ型

Adjustable Ankle Orthosis



## Indication:

- Use for fixation of foot and ankle fracture, lower tibiaand fibula fracture.

## Characteristic:

- Adjustable button is easy to operate which can be locked in a certain angle or limited range and can be adjusted protected range-of-motion according to patient's condition of rehabilitation;

- ☒•Anti-slip design of sole, enhance the stability

- ☒• Universal design fits right and left foot

## Size:

Size	Maximum foot length(cm)	Shoe size
S	24.5	36-40
M	26.5	40-42
L	29	42-45



TJ-603 (1)

跟腱靴固定器 I 型

Short Knee Orthosis



## Indication:

- Use for foot and ankle ligament, soft tissue and tendon injury conservative treatment.
- After Forefoot ,toes and ankle's injury conservative treatment.
- Conservative treatment of after operation of distal fibula fracture , stead forefoot, toes fracture .

## Characteristic:

- The lining is made of skin-friendly fabric, which is soft and comfortable;
- Double side airline protection, fit to protect Achilles tendon;
- The push-type airline is easy to use and operate.

## Size:

Size	Maximum foot length(cm)	Shoe size
S	24.5	36-40
M	26.5	40-42
L	29	42-45

TJ-603 (2)

跟腱靴固定器 V 型

Short Knee Orthosis



## Indication:

- Use for foot and ankle ligament, soft tissue and tendon injury conservative treatment.
- After Forefoot ,toes and ankle's injury conservative treatment.
- Conservative treatment of after operation of distal fibula fracture , stead forefoot, toes fracture .

## Characteristic:

- The lining is made of skin-friendly fabric, which is soft and comfortable;
- Double side airline protection, fit to protect Achilles tendon;
- The push-type airline is easy to use and operate.

## Size:

Size	Maximum foot length(cm)	Shoe size
S	24.5	36-40
M	26.5	40-42
L	29	42-45

TJ-603 (3)

跟腱靴固定器VII型

Short Knee Orthosis



## Indication:

- Use for foot and ankle ligament, soft tissue and tendon injury conservative treatment.
- After Forefoot ,toes and ankle's injury conservative treatment.
- Conservative treatment of after operation of distal fibula fracture , stead forefoot, toes fracture .

## Characteristic:

- The lining is made of skin-friendly fabric, which is soft and comfortable;
- Double side airline protection, fit to protect Achilles tendon;
- The push-type airline is easy to use and operate.

## Size:

Size	Maximum foot length(cm)	Shoe size
S	24.5	36-40
M	26.5	40-42
L	29	42-45

TJ-604

小腿超踝固定带

Ultra Ankle Fixation Brace



### **Indication:**

- Tibia and fibula fracture fixation, ankle fracture fixation, tibial fibular ligament injury.

### **Characteristic:**

- The brace overall wraps leg and foot, which lined with aluminum alloy support plate to ensure fixed effect.

### **Size:**

Universal size



TJ-605

踝足固定器 I 型

Night Splint



## Indications:

- Used for tendonitis,plantar fasciitis, pain caused by ankle and foot strain.
- Prevention of foot drop and Achilles tendon contracture for bedridden patients;
- Fixation for ankle fracture, fracture of lower tibia and fibula and ankle ligament injury.

## Characteristic:

- The product avoids secondary damage.The soft pad ensures maximum comfort on fixing achilles tendon.

## Size:

Size	Foot Length(cm)
S	<20
M	21-27
L	>28

TJ-606

踝关节固定带 I 型

Ankle Brace



## Indications:

- Conservative treatment of ankle joint and soft tissue injury.

## Characteristic:

- Two external straps wrap to enhance stability.

## Size:

Universal size

# TJ-606 (1)

## 踝关节固定带 II 型

### Ankle Brace



### Indications:

- Conservative treatment of the ankle joint and its soft tissues,ankle ligament injuries, ankle instability, etc.

### Characteristic:

- Knob type fine-tuning locking system, easy to wear, uniform pressure, easy adjustment of tightness;
- Plastic support strips conform to physiological design, tightly fit, and secure;
- High-quality fabrics, skin-friendly lining, moisture wicking, comfortable and breathable;
- Excellent workmanship, lightweight and short, stylish appearance.

### Size:

Size
S (Left/Right)
L (Left/Right)

TJ-606 (2)

踝关节固定带 V 型

Ankle Brace



## Indications:

- Conservative treatment of the ankle joint and its soft tissues,ankle ligament injuries, ankle instability, etc.

## Characteristic:

- Plastic support plates on both sides are designed according to the ankle curve to provide support for the ankle.
- Using drawstring design, easy to wear and adjust .

## Size:

Size	Foot length(cm)	Shoe size
S	22.5-24	36-40
M	24.5-26	40-42
L	26.5-28.5	42-45
XL	>28.5	>47



## TJ-606 (3) 踝关节固定带VI型 Ankle Brace



### Indications:

- Ligament injury, acute ankle injury, postoperative rehabilitation

### Characteristic:

- Plastic support plates on both sides are designed according to the ankle curve to provide support for the ankle.
- Shoelace design, easy to adjust

### Size:

Size	Foot length(cm)	Shoe size
S	22.5-24	36-40
M	24.5-26	40-42
L	26.5-28.5	42-45
XL	>28.5	>47

TJ-607 (1)

踝关节固定带IV型

Ankle Splint



## Indications:

- Chronic instability of the ankle joint, varus and valgus sprain, fibular ligament injury, postoperative care and fixation of the ligament, promote the absorption of hematoma and edema, and walk protection during rehabilitation.

## Characteristic:

- Ergonomic support board, more fit, comfortable and stable.
- Adjustable two-way circuitous strap, convenient adjustment, easy to wear, hard to loose.
- Adjustable heel straps can freely adjust the tightness, making it suitable for different sizes of ankles.
- 3 sets of pads can option (airbag, foam, gel).

## Size:

S、L

TJ-607 (2)

踝关节固定带Ⅷ型

Ankle Brace



## Indications:

- Suitable for ankle sprain, ankle ligament injury, habitual sprain, sports crowd protection

## Characteristic:

- The material is light and breathable to make it more comfortable to wear.
- Ergonomic design, provide solid support, reduce the risk of injury.
- The single strap design is convenient to wear and the cross compression fixing is more reliable.
- It fits the ankle height and can wear shoes normally at the same time to help rehabilitation exercise.

## Size:

Size	Foot length(cm)	Shoe size
S(Left/Right)	22.5-24	35-38
M(Left/Right)	24.5-26	39-42
L(Left/Right)	26.5-28.5	43-47

TJ-607 (3)

踝关节固定带VII型

Ankle Brace



## Indications:

- Ligament injury, acute ankle injury, postoperative rehabilitation.

## Characteristic:

- The foot strap can be automatically adjusted according to the width of the foot, and to stabilize the ankle without displacement.
- Through the plantar on the instep of the figure-eight winding fixed, stabilize the ankle and strengthen protection.
- The foam filled air pad wraps and squeezes both sides of the ankle muscles gently, and relieve discomfort during exercise.
- The shell fully supports and protects the ankle.

## Size:

Size	Foot length(cm)	Shoe size
S(Left/Right)	22.5-24	35-38
M(Left/Right)	24.5-26	39-42
L(Left/Right)	26.5-28.5	43-47
XL(Left/Right)	>28.5	>47



TJ-701

卷式夹板

Roll Splint



## Indications:

- Superior first aid product for fractured or injured limbs,which is also the necessary need for variety of first aid kits.

## Characteristic:

- It is made of a thin core of metal plate, sandwiched between two layers of foam.
- Small in size, waterproof and easy to use.
- Achieve the best support simply shaped to any curves by hand.

## Size:

Length:20-100cm

Width:10-15cm